

RCC-ECOSTAND-02R

Height Adjustable EcoStand compatible with Power on Pole System(HA)



Product Features

- Easy cleaning
- Adjustable or Fixed height
- Small footprint for space saving
- Infection control friendly
- Lightweight & maneuverable
- Compatible with our GeniTec™ power on pole system
- Monitor and device mounts available

Specification

SPEC

Height	1185 mm to 1585 mm - Adjustable 400mm
Max Tilt Angle	10°(in a static tilt test) As per EN 60601-1 section 9.4.2.1
Base size	550mm
Accessories Available	Handle, Basket, Medi-Rail, Cable Hook, Drip Pole, PSU Adapter Holder, Keyboard Tray
Mounting Options	Single and double accessory pole mounts available (For attaching accessories)
Geni-Tec System	P.O.P (Power-On-Pole) Single & Dual Dock available
Castors	5 castors, 75mm, All braked
Compliance	CE, UKCA

Ordering Information

Ordering Information

RCC-ECOSTAND-02R (HA), Height Adjustable EcoStand compatible with Power on Pole System(400mm)/Extrusion length (1000mm)/5 Star base 600mm/Standard 75mm Wheels

RCC-ECOSTAND-01R (FH), Fixed Height EcoStand compatible with Power on Pole System/Extrusion length (1000mm)/5 Star base 600mm/Standard 75mm Wheels