

RCC-MEDISTAND-01R

Fixed Height Medi-Stand (FH)



Product Features

- Seamless moulding and smooth surfaces for easy cleaning
- Height adjustable Fixed height
- Power ready with our GeniTec™ power system for power on the go
- Lightweight maneuverability
- Integrated cable management
- Small footprint for space saving

Specification

SPEC

Height	Various (Standard length 1250mm)
Max Tilt Angle	10 Degrees (in a static tilt test) As per EN-60601-1 section 9.4.2.1
Base size	675mm Diameter - 5 wheel base 561mm x 557 mm - 4 wheel base
Accessories Available	Handle, Basket, Medi-Rail, Cable Hook, Drip Pole, PSU Adapter Holder, Keyboard Tray
Castors	5 castors, 100 mm, all braked 4 castors, 100 mm, all braked
Max Mount Weight	12 kg (Alternatives Available)
Max Carrying Load	18 kg (subject to height)
Compliance	CE / UKCA

Ordering Information

Ordering Information	RCC-MEDISTAND-01R (FH), Fixed Height Medi-Stand, Height various (Standard length 1250mm)
	RCC-MEDISTAND-02R (HA), Height Adjustable Medi-Stand, Height 1250mm-1675mm, range: (400mm Adjustable) Or (300mm optional)