

RCC-GENISTAND-01R

Height Adjustable Geni-Stand with Single or Dual Power Management System (HA)



Product Features

- Easy Cleaning
- Ergonomically Designed
- Heavy Duty Castors, all braking
- Custom mount fittings and accessories
- Height adjustable option available
- Lightweight for easy maneuverability
- Small footprint for space saving
- Robust construction for busy environments
- Cable management system, no trip hazards
- 216 Wh power module for extended operation between charges

Specification

SPEC

Height	1155mm to 1555mm - Adjustable 400mm
Max Tilt Angle	10°(in a static tilt test) As per EN 60601-1 section 9.4.2.1
Base size	675mm Diameter
Accessories Available	Handle, Basket, Medi-Rail, Cable Hook, Drip Pole, PSU Adapter Holder, Keyboard Tray
Mounting Options	Accessory adaptor mounts compatible (channel adaptor)
Geni-Tec System	Single & Dual-dock systems available
Castors	5 Castors, 100mm, All braked 4 Castors, 100mm, All braked
Compliance	Stand: CE, UKCA Power System: EN/ES/IEC 60601-1, CE, FCC Power Module: IEC 62133, UL 2054, UN38.3